

Operation Manual

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2.3 QC Material LOT

2.3.1 QC Material LOT	Depending on the produced QC Material LOT, it consists of Made Date, Serial, Expire Date, etc. An ITEM file is required to create a QC Material LOT.
2.3.2 QC Material LOT file	QC Material LOT XML file that stores LOT information and a Barcode PDF file are created. Ex) VEUDx_QC_LOT_TBI-assay_230921001.xml VEUDx_QC_LOT_TBI_230921001_BarCode.pdf

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

- You must input information about the name, unit, and RSMP Length of the markers to be used for the item.
- Based on the input information, the fluorescence image is detected by RSMP Length and the result is calculated.

VEUDx Experiment Editor V1.1.31.8 - Make ITEM

ITEM Summary Protocol Home Setting

ITEM Name TBI-assay Initialize

Length	Marker	Use ExtMFI	Count	Unit	QC Conc.			LOD	minMFI	Cal MFI			
					CutOff	QC1	QC2			QC1 min	QC1 Max	QC2 min	QC2 Max
130		<input type="checkbox"/>		pg/ml									
200	UCH-L1	<input type="checkbox"/>		pg/ml	138	200	50	0.78	1500	5000		3000	
250		<input type="checkbox"/>		pg/ml									
300	GFAP	<input type="checkbox"/>		pg/ml	24	100	25	0.39	1500	5000		3000	
350		<input type="checkbox"/>		pg/ml									
400		<input type="checkbox"/>		pg/ml									
450		<input type="checkbox"/>		pg/ml									
500		<input type="checkbox"/>		pg/ml									

Total Estimated Time (Min.) 58 QC Replication Not Specified

+ RT 6 samples assumptions Fluorescence Exposure Time1(ms) 500

Reacting Time Only 49 Fluorescence Exposure Time2(ms) 0

Pixel Cut (%) Dilution Factor 4 ☐ Use Extended MFI

Bottom 25 ☐ Make RSMP QC Data

Top 5 Buffer Well None ☐ Optical photography only

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'Fluorescence Exposure' value.

(High MFI values can affect fluorescence saturation)

- Optical photography only

Check if you only want to image capture and Detecting

(RSMP is demagnetized)

- Optical + Demagnetization only

Check if you only want to image capture and Detecting


(RSMP is not demagnetized)

(RSMP should be in imaging well 8)

- Make RSMP QC Data

If you select the 'Make RSMP QC Data' Option, you can acquire MFI, CV(%) values for each RSMP and for the entire RSMP.

(Used to check RSMP Ab coupling QC status)

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
- Cartridge Barcode : VEUDx-C/V3/230922001

- Reaction Protocol
 1 : Heating | 39°
 2 : Reacting | Well 2 | Shuffle every 5 min | 30 min
 3 : Washing | Well 3 | 30 Secs x 2
 4 : Reacting | Well 4 | Shuffle every 2 min | 5 min
 5 : Washing | Well 6 | 30 Secs x 2
 6 : Reacting | Well 5 | Shuffle every 2 min | 5 min
 7 : T-Tip Cleaning | Well 1 | 30 secs
 8 : Washing | Well 7 | 20 secs x 2

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
6.1 Select ITEM file


6.1.1 Select ITEM file

Select ITEM file to make QC Material LOT.

6.2 Make QC Material LOT

6.2.1 Start

The ITEM name read from the ITEM file are displayed.

- 6.2.2 Enter QC Material LOT information**
- QC Material LOT creation date
 - QC Material LOT Serial
 - QC Material LOT Expire Date can be entered.
- *Barcode Name, Note can be used as a reference.

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
6.2.4 Save QC Material LOT Click the "Save QC LOT(+PDF)" button
Two QC LOT XML files and Barcode PDF files are created.

- QC Material LOT XML file (ex, VEUDx_QC_LOT_TBI-assay_230922001.xml)
used to save work

- QC Material Barcode PDF file
(ex, VEUDx_QC_LOT_TBI-assay_230922001_BarCode.pdf)

It is used for reading VEUDx equipment QC Material Barcode by printing it out on paper.

< VEUDx equipment QC Material LOT reading screen >

< PDF file for QC Material LOT Barcode output >

6.2.5 Open QC Material LOT Editing is possible by selecting the saved LOT XML by pressing the "Open QC LOT" button.

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